

# AL001

## UNIVERSAL MANUAL CARRIER



### AFFORDABLE, FULLY-MANUAL SOLUTION

Economical, fully manual solution carries virtually any power wheelchair or scooter.

## AL001 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	500 lb
<b>PLATFORM SIZE</b>	28" x 48" platform with a 60" folding ramp
<b>HITCH CLASS</b>	II or III
<b>AVAILABLE OPTION(S)</b>	- Cover - Tie-down straps

<b>INSTALLED WEIGHT</b>	84 lb
<b>ACCOMMODATES</b>	- Manual wheelchairs - Power wheelchairs - Scooters
<b>WARRANTY</b>	1-year



Harmar.com  
800.833.0478

08AUG2020 | MKT-000129 REV B

